3diag - U-B2m - T1A

REF TD-42531-U

LOT **531U-10**

2025-01

UDI

(01) 08434477302060
(10) 531U-10
(17) 250130
(240) C33726

EN Beta-2 Microglobulin - Urine - for Turbidimetry

INSTRUCTIONS FOR USE

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IT Beta-2 Microglobulina - Urina - per Turbidimetria

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<https://www.3diag.com/002>

for IFU, Scan or follow link,
and select LOT

TRIMERO Diagnostics, SL

c. València 558, 4t 2a - 08026 Barcelona (Spain)
☎ +34 93 244 86 79 - www.3diag.com

INSTRUCTIONS FOR USE

Reagents for professional use,
for *In Vitro* use only in clinical laboratory (IVD)

3diag - U-B2m - TIA

Beta-2 Microglobulin - Urine
for Turbidimetry

REF TD-42531-U

INTENDED USE

Quantitative determination of Beta-2 Microglobulin (B2m), in human urine, by turbidimetric method in automatic Clinical Chemistry Analyzers.

PRINCIPLE OF THE METHOD

The specific antibodies (Ab) of the reagent, bound to polystyrene particles, when combined with the antigens (Ag) of the patient sample, form insoluble compounds causing a change in the absorbance and dispersion of the light, proportional to the antigen concentration, which can be quantified by turbidimetric (TIA) or nephelometric (NIA) method, by comparison with calibrators of known concentration.

CONTENTS - COMPOSITION - PREPARATION

- Antiserum Reagent: **REAG** **Ab** **U-B2m**
REF TD-42531-RA ▽ 100 test^(*) - 3 ml
Anti-human B2m antibodies bound to polystyrene particles.
- Reaction Buffer: **BUF** **U-B2m**
REF TD-42531-BF ▽ 100 test^(*) - 15 ml
TRIS Buffer, with PEG.

Note (*1): with the recommended general assay parameters.

The reagents are ready for use and require no preparation.

Before each use it is convenient that the reagents are homogenized, shaking them gently avoiding the formation of foam or bubbles.

As a preservative, the reagents contain <0.1% (1 g/l) Sodium Azide (NaN₃).

WARNINGS - PRECAUTIONS

- Sodium Azide is toxic. Even if sodium azide is not harmful at the concentration present in the reagents, take the necessary precautions to avoid accidental ingestion or contact with the eyes.
- Sodium Azide can react with lead or copper to give explosive compounds. For disposal it is recommended to rinse with plenty of running water to avoid accumulation in drains.
- Since the absence of infectious agents can not be proven with absolute certainty, components containing materials of human or animal origin must be handled with caution, as potentially infectious, following the recommended safety standards for biological risk.
- Do not mix components belonging to different lot kits.
- Clinical diagnosis should not be based on the results of a single test, but should always integrate all relevant clinical and laboratory data.

STORAGE - SHELF LIFE

- Store refrigerated at +2...+8°C. Do not freeze, as the functionality of the reagents may be altered.
- Properly stored and unopened, the reagents are stable until the expiration date indicated on the label.
- Once opened, the shelf life of the reagents is at least 4 weeks, provided that after each use they are stored immediately in the original containers, tightly capped and refrigerated at +2...+8°C. This information should be taken as a guideline given that, obviously, the shelf life depends on the particular environmental and use conditions, which may differ from those of the stability studies carried out.

MATERIALS NEEDED, NOT SUPPLIED

- Automatic Clinical Chemistry Analyzer, capable of running photometric assays at 540...660 nm, and accessories: reagent containers, cuvettes, etc..
- | | |
|------------------------------|---------------------|
| • 3diag - U-B2m - CAL SET | REF TD-42532 |
| • 3diag - B2m - CONTROL (x3) | REF TD-42523 |
| • 3diag - B2m - CONTROL | REF TD-42533 |

SAMPLES

Fresh urine.

It is usual the use of a 24-hour urine aliquot, although in the literature some authors recommend the use of random urine, preferably the second morning void, and expressing the results relative to urinary creatinine.

Prior to the analysis, the samples should be centrifuged until a clear and transparent supernatant is obtained⁽¹⁾.

For the determination of specific proteins, centrifugation of urine samples at 3000⁽²⁾-5000⁽³⁾ g for 10 minutes is the standard practice in the laboratory.

In bibliography⁽⁴⁾ it is reported a stability in urine of:

- Refrigerated: 48 hours
- Frozen: 60 days.

The pH must be adjusted between 6 and 8, since the B2M degrades rapidly if pH <6. It is recommended⁽⁴⁾ to use 1M NaOH.

PROCEDURE

If necessary, carefully transfer the reagents to the containers used by the analyzer, preventing leakage and foaming or bubbles.

To program and calibrate assays, follow the instructions for use of the analyzer used, with the recommended general parameters that are detailed below. Please, contact the Customer Support Service (✉ support@3diag.com - ☎ +34 93 244 86 79) for further information about applications to specific analyzers.

Assay Parameters

- ①Dispense and mix:
 - Sample/Calibrator/Control: 12 µl (diluted 1:5)
 - **BUF** **U-B2m** 150 µl
- ②Incubate a fixed time between 1 and 5 minutes
- ③Dispense and mix:
 - **REAG** **Ab** **U-B2m** 30 µl
- ④Read absorbance A1 (Blank) at 540...660 nm
- ⑤Incubate a fixed time of between 3 and 5 minutes
- ⑥Read absorbance A2 (Final Point) at 540...660 nm
- ⑦Interpolate the absorbance increment (A2-A1) of the samples and controls in the curve obtained with the calibrators
- ⑧Samples with concentrations higher than the upper limit of the assay range should be analyzed again, diluted manually with Physiological Solution, or by programming a larger sample dilution in the analyzer, to recover a value close to the midpoint of the measurement range.

As an alternative, reagents can be mixed as first step, and the sample dispensed as starter.

Calibration Parameters

- Use the 3diag - U-B2m - CAL SET.
- If the analyzer allows it, it is recommended to program two replicates of each calibration point.

- The calibration is Non-linear. For the calculation it is recommended to use a 3rd Order Polynomial, a Logit or a Polygonal adjustment.
- The assay must be recalibrated, at least when a new batch of reagents is used or when its parameterization is changed.

PERFORMANCES OF THE METHOD

Detailed information on the characteristics and performances of the assay is given in the Technical Report, available on the website (www.3diag.com) or upon request to the Customer Support Service (✉ support@3diag.com - ☎ +34 93 244 86 79).

QUALITY CONTROL

To monitor performances, it is recommended that internal controls be inserted into each analytical series. It is recommended to use the controls of **3diag - B2m - CONTROL**.

Each laboratory should establish its own quality scheme and corrective actions if controls do not meet the assigned tolerances.

The reagents have been subjected to quality control checks and should react as described in these instructions. Therefore, as a general recommendation, in case the controls do not give the expected reaction, as a precaution all reagents should be considered unreliable until their operation has been checked.

TRACEABILITY

Values are referred to the *1st International Standard for Beta2 Microglobulin* (code: *B2M*) of the *WHO (World Health Organisation)*.

REFERENCE INTERVALS

It is always advisable for each laboratory to establish its own reference values.

The bibliography reports that a B2m excretion greather than 0.37 mg/24h⁽⁵⁾ or a concentration of 0.3 mg/l⁽⁴⁾ can be considered significant.

CLINICAL SIGNIFICANCE

Elevated serum levels of B2m are either due to an increased production or to a decrease in the glomerular filtration rate which determines its accumulation. The increased rate of synthesis may be caused by many different diseases, in particular those of the lymphoid immune system, hematologic malignancies, viral infections, autoimmune diseases and others. B2m is hence not used as a specific diagnostic marker for any pathology in particular, but is rather widely used as a quantitative prognostic marker in many diseases. To take just one example, measurement of B2m levels is included in the International Staging System for Multiple Myeloma.

Measurement of B2m concentrations in urine is a very sensitive marker of tubular dysfunction due to any cause or pathology. Although as a tubular marker it suffers from the problems of instability in acidic environments (pH<6) and variable production levels.

SYMBOLS

In addition to the harmonized symbols provided on the European Standard EN 980:2008, in the labels and instructions of use has been used the complementary symbology proposed⁽⁶⁾ by the *EDMA (European Diagnostic Manufacturers Association)*, whose meaning is detailed below.

REAG	Reagent
Ab	Antibody / Antiserum
BUF	Buffer
U-B2m	Beta-2 Microglobulin - Urine

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I valori sono riferiti al *1st International Standard for Beta2 Microglobulin* (code: *B2M*) della *WHO (World Health Organisation)*.

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
NOTE #1: Applications available for **Optilite®**, **Alinity c**, **Architect c**, **ADVIA®** series, **AU®** series, **cobas®** series, and other analyzers by request

IMAGE® and **AU®** are registered trademarks of Beckman Coulter, Inc, Fullerton, CA.
BN™ and **ADVIA®** are registered trademarks of Siemens Healthcare Diagnostics Products GmbH, Marburg, Deutschland.
Optilite® is a registered trademark of The Binding Site Group Ltd, Birmingham, U.K.
Alinity and **Architect** related brand marks are registered trademarks of Abbott Laboratories, Abbott Park, Illinois, USA
cobas® and related brand marks are trademarks of Roche Diagnostics Ltd, Rottkreuz, Switzerland

ALL PROTEINS INCLUDED IN  **Accuracy 365** Online Quality Control Further information and registration at www.accuracy365.com

Accuracy 365 aims to monitor and compare, between laboratories, the results obtained in the Internal Quality Control processes of end-users of specific protein control materials manufactured by **TRIMERO Diagnostics**.

Main features:

- 100 % free and lifetime.
- Cloud application, Available 24 hours, 365 days per year.
- Intuitive interface and very simple to use.
- Calculation and representation of statistical data in real time.
- Anonymous participation, data are only accessible by its owner.

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- (2) "Alpha-1-Microglobulin (A1M) - IMMAGE® Immunochemistry Systems Chemistry Information Sheet", © Copyright 2017 Beckman Coulter, Inc..