

## 3diag - C1In - 800

**C1 (Esterase) Inhibitor - for *IMMAGE*<sup>®</sup> 800**

**REF** TD-42581

**IVD** CE

## ANNEX to Instruction for Use:

**Application for *IMMAGE*<sup>®</sup> and *IMMAGE*<sup>®</sup> 800 Immunochemical Systems**

<b><u>PROTOCOL INFORMATION</u></b>			
Reagent Lot Number:	(*1)	Units:	mg/dL or mg/L
Cartridge Lot Number:	(*1)	Protocol:	Non-competitive Nephelometry
Cartridge Serial Number:	(*1)	AGXS Limit:	not defined
Reagent Expiry date:	(*1)		
Sample or Dilution Volume:	20 µL	Diluent:	DIL1 (*2)
Reaction Buffer Volume:	195 µL	Buffer:	BUF1 (*3)
Calibration Dilution:	1 : 7	Gain:	3
Compartment A Volume:	20 µL	Sample Dilution:	1 : 7
Compartment B Volume:	0 µL (not used)	Reaction Time:	2.5 minutes
<b><u>CALIBRATION INFORMATION</u></b>			
Levels:	6	Replicates:	2
Update Level:	4 (*4)	Replicates:	2 (*4)
Model:	Logistic 4 parameters or 3rd order Polynomial (select the best adjustment)		
<b><u>NOTES</u></b>			
<p>(*1): See the labels of the Reagent Cartridge.</p> <p>(*2): Beckman Coulter P/N no. 447640. For the <i>IMMAGE</i><sup>®</sup> it has to be defined as DIL10.</p> <p>(*3): Beckman Coulter P/N no. 447650. For the <i>IMMAGE</i><sup>®</sup> it has to be defined as BUF10.</p> <p>(*4): The use of the calibration update, with only one calibrator level, is discouraged.</p>			