

## N Latex C1q Kit

### ANNEX to IFU: *Atellica® NEPH 630 / BN ProSpec® Systems* - Proposal of Application

#### REF TD-42550 - C1q Complement - for *BN™ Series and Atellica® NEPH 630*

*BN™, ProSpec® and Atellica®* are registered trademarks of *Siemens Healthcare Diagnostics Products GmbH, Marburg, Deutschland*

<u>Assay-Definition</u>					
<b>Assay-ID:</b>	200	<b>Name:</b>	C1q_TD	<b>long name:</b>	C1q
<b>Version:</b>	10721	<b>Assayclass:</b>	Neph	<b>Assaygroup:</b>	Serum
<b>Result unit:</b>	mg/dl	<b>lower reference range</b>	0	<b>upper reference range</b>	9999999999
<b>Sample dilution:</b>	1:100	<b>Min. dilution:</b>	1:20		
<b>Max. dilution:</b>	unassigned				
 <b><u>Evaluation:</u></b>					
<b>Measurement duration [1/10 sec.]:</b>	7200	<b>Measurement transfer:</b>		2	
<b>Evaluation method:</b>	9 - Two point (BNP)				
<b>Measurement window [1/10 sec.] from:</b>	100	<b>until:</b>	7200		
<b>Blockmode:</b>	0	<b>Overflow:</b>	0		
 <b><u>Calibration:</u></b>					
<b>Calibration method:</b>	LogitLog	<b>Deviation [%]:</b>	10	<b>Calibration valid until:</b>	0
<b>Calibrator:</b>	6705 - C1Q CAL L				
<b>Start dilution:</b>	1:10	<b>Supporting points:</b>	5		
 <b><u>Special parameters:</u></b>					
<b>Parameter 1:</b>	2	<b>Parameter 2:</b>	0,5	<b>Parameter 3:</b>	-30
				<b>Parameter 4:</b>	80
 <b><u>Controls</u></b>					
<b>Control-ID</b>	6706	<b>Name</b>	C1q Control H		
<b>Dilution</b>	1:100	<b>Deviation [%]</b>	20	<b>Level</b>	High
<b>Control-ID</b>	6707	<b>Name</b>	C1q Control L		
<b>Dilution</b>	1:100	<b>Deviation [%]</b>	20	<b>Level</b>	Low
 <b><u>Transfers:</u></b>					
<b>Transfer number:</b>	1				
<b>Mixing:</b>	normal	<b>Start wash:</b>	13 - intensive	<b>End wash:</b>	06 - goto waste
<b>Incubation time [1/10 sec.]:</b>	0	<b>Incubation delta [1/10 sec.]:</b>	0		
 <b><u>Medium:</u></b>					
<b>Medium number:</b>	1				
	<b>Type:</b>	SystemLiquid	<b>Medium ID:</b>	Buffer	
	<b>Syringe profile:</b>	3	<b>Volume [µl]:</b>	165	
<b>Medium number:</b>	2				
	<b>Type:</b>	SystemLiquid	<b>Medium ID:</b>	Air bubble	
	<b>Syringe profile:</b>	2	<b>Volume [µl]:</b>	2	
<b>Medium number:</b>	3				
	<b>Type:</b>	Reagent	<b>Medium ID:</b>	6796 - C1q Enh	
	<b>Syringe profile:</b>	3	<b>Volume [µl]:</b>	40	
<b>Medium number:</b>	4				
	<b>Type:</b>	Sample	<b>Volume [µl]:</b>	20	
<b>Medium number:</b>	5				
	<b>Type:</b>	Destination	<b>Volume [µl]:</b>	-127	
	<b>Syringe profile:</b>	3			

## N Latex C1q Kit

**ANNEX to IFU: *Atellica® NEPH 630 / BN ProSpec® Systems* - Proposal of Application**

Assay-Definition

<b>Assay-ID:</b>	200	<b>Name:</b>	C1q_TD	<b>long name:</b>	C1q
<hr/>					
<b>Transfer number:</b>	2				
<b>Mixing:</b>	normal	<b>Start wash:</b>	13 - intensive	<b>End wash:</b>	06 - goto waste
<b>Incubation time [1/10 sec.]:</b>			0	<b>Incubation delta [1/10 sec.]:</b>	0
<u>Medium:</u>					
<b>Medium number:</b>	1	<b>Type:</b>	SystemLiquid	<b>Medium ID:</b>	Diluens
		<b>Syringe profile:</b>	7	<b>Volume [µl]:</b>	95
<b>Medium number:</b>	2	<b>Type:</b>	SystemLiquid	<b>Medium ID:</b>	Air bubble
		<b>Syringe profile:</b>	2	<b>Volume [µl]:</b>	2
<b>Medium number:</b>	3	<b>Type:</b>	Reagent	<b>Medium ID:</b>	6703 - C1q Ab
		<b>Syringe profile:</b>	2	<b>Volume [µl]:</b>	20
<b>Medium number:</b>	4	<b>Type:</b>	Destination		
		<b>Syringe profile:</b>	6	<b>Volume [µl]:</b>	-117