

## 3diag - C1q - TIA

### ANNEX to IFU: *Optilite*<sup>®</sup> - AOR Test Parameters - Proposal of Application

#### **REF** TD-42551 - C1q Complement - For Turbidimetry

*Optilite*<sup>®</sup> is a registered trademark of *The Binding Site Group Ltd.*, Birmingham, U.K.

#### INFO

<i>Name</i>	User Defined		
<i>Tag</i>	Read-only field		
<i>Version number</i>	Updated when new settings are saved		
<i>Full name</i>	User Defined		
<i>Online name</i>	User Defined		
<i>Type</i>	<b>Photometric</b>	<i>Number of decimals</i>	<b>2</b> (minimum recommended 1)
<i>In use</i>	<b>Yes</b>	<i>Correction factor</i>	<b>1.000</b>
<i>Acceptance</i>	<b>Manual</b>	<i>Correction bias</i>	<b>0.000</b>
<i>Result Unit</i>	<b>mg/dl</b>		
<i>Sample Type</i>	<b>Serum</b>		
<i>Test version ID</i>	Read-only field		
<i>Last time changed</i>	Date and time of last changes		
<i>User name</i>	User ID, who modified the definition		

#### FLOW

<i>Blank type</i>	<b>Yes</b>	<i>Primary dilution 1+</i>	<b>29</b> (Read-only field)
		<i>Dispensed volume</i>	Read-only field

#### 1st Step - Reagent

<i>Reagent</i>	Select Reagent from the drop down menu - Reaction Buffer - Use <b>REF</b> TD-42551-BF - <b>BUF C1q</b>		
<i>Volume (µl)</i>	<b>120</b>		
<i>Dispense with</i>	<b>Extra</b>		
<i>Extra volume (µl)</i>	<b>10</b>		
<i>Syringe speed</i>	<b>Medium</b>		
<i>Replacing Reagent</i>	<b>None</b>		

#### 2nd Step - Reagent

<i>Reagent</i>	Select Reagent from the drop down menu - Antiserum Reagent - Use <b>REF</b> TD-42551-RA - <b>REAG Ab C1q</b>		
<i>Volume (µl)</i>	<b>30</b>		
<i>Dispense with</i>	<b>Extra</b>		
<i>Extra volume (µl)</i>	<b>10</b>		
<i>Syringe speed</i>	<b>Slow</b>		
<i>Replacing Reagent</i>	<b>None</b>		

#### 3rd Step - Mix

#### 4th Step - Incubation

<i>Time (sec)</i>	<b>54</b>	<i>Actual time (sec)</i>	Read-only field
-------------------	-----------	--------------------------	-----------------

## 3diag - C1q - TIA

### ANNEX to IFU: *Optilite*<sup>®</sup> - AOR Test Parameters - Proposal of Application

<b>5th Step - Sample</b>				
Volume (µl)	<b>6</b>			
Dispense with	<b>Extra</b>			
Extra volume (µl)	<b>10</b>			
Extra wash	<b>No</b>			
<b>6th Step - Mix</b>				
<b>7th Step - Incubation</b>				
Time (sec)	<b>18</b>	Actual time (sec)		Read-only field
<b>8th Step - End-point (Blank)</b>				
Blank resp. min (A)	<b>0</b> (Not used)	Blank resp. max (A)		* (Not used)
<b>9th Step - Incubation</b>				
Time (sec)	<b>297</b>	Actual time (sec)		Read-only field
<b>10th Step - End-point (Measurement)</b>				
Main wavelength	<b>600</b>	Side wavelength		<b>None</b>
Residual net abs. (A)	<b>0</b> (Not used)	Delta abs. check min. (A)		* (Not used)
<b><u>DILUTION</u></b>				
Dilution with	<b>Diluent</b>			
Primary dilution 1+	<b>29</b>			
<b>Neat sample</b>				
Dispense with	<b>Extra</b>			
Volume (µl)	<b>10</b>			
<b>Diluent</b>				
Sample diluent ID	<b>Diluent 1</b>			
Calibrator diluent ID	<b>Diluent 1</b>			
Dispense with	<b>Extra</b>			
Volume (µl)	<b>10</b>			
<b><u>LIMITS</u> ( User defined, the proposed parameters have only value as a recommendation )</b>				
	<b>Measuring range</b>		<b>Next dilution ratio 1+</b>	
	<b>Min</b>	<b>Max</b>	<b>Low</b>	<b>High</b>
Primary dilution	<b>See note (*1)</b>	<b>See note (*2)</b>	<b>9</b>	<b>119</b>
2nd, 3rd & 4th dilution	* (Not used)	* (Not used)	* (Not used)	* (Not used)
<b>Test limit</b>				
Critical limit	<b>See note (*3)</b>	<b>See note (*4)</b>		
Init. abs.	* (Not used)	* (Not used)		
Reference ranges	<b>User Defined</b>			

## 3diag - C1q - TIA

### ANNEX to IFU: *Optilite*<sup>®</sup> - AOR Test Parameters - Proposal of Application

#### CALIBRATION

Calibration type	<b>Logit-Log4</b>	Abs. error (A)	* (Not used)
Repeat time (days)	<b>0</b> (Not used, see note (*5))	Rel. error (%)	* (Not used)
Points / calibrator	<b>Duplicate</b>	Factor limit min.	* (Not used)
Acceptance	<b>Manual</b>	Factor limit max.	* (Not used)
Concentration axis	<b>Linear</b>	Bias limit min.	* (Not used)
Response axis	<b>Linear</b>	Bias limit max.	* (Not used)

Calibrators ( User-defined Calibrators are defined in "F4 > Cal/Ctrl definition" - as Calibrator, use REF TD-42552 - **3diag - C1q - CAL** )

Nbr	Calibrator	Current Lot	Concentration	Dilution 1+
1	Select Cal.	Read-only value	Read-only value	<b>199</b>
2	Select Cal.	Read-only value	Read-only value	<b>119</b>
3	Select Cal.	Read-only value	Read-only value	<b>59</b>
4	Select Cal.	Read-only value	Read-only value	<b>44</b>
5	Select Cal.	Read-only value	Read-only value	<b>29</b>

#### NOTES

- (\*1) Equal to lowest calibration point (Nbr 1) concentration  
 ( Calculated as:  $\text{Cal. Concentration} * (\text{Sample Dilution} / \text{Cal Nbr 1 Dilution}) = \text{Cal. Concentration} * (30 / 200)$  )  
 ( If **AUTOM** flag is selected, limits are defined and recalculated from the calibration )
- (\*2) Equal to highest calibration point (Nbr 5) concentration  
 ( Calculated as:  $\text{Cal. Concentration} * (\text{Sample Dilution} / \text{Cal Nbr 5 Dilution}) = \text{Cal. Concentration} * (30 / 30)$  )  
 ( If **AUTOM** flag is selected, limits are defined and recalculated from the calibration )
- (\*3) Equal to minimum measuring range divided by the re-concentration factor of the primary dilution  
 ( Calculated as:  $\text{Min. Meas. Range} / (\text{Sample Dilution} / \text{Low Next Dil. Ratio}) = \text{Min. Meas. Range} / (30 / 10) =$   
 $= \text{Cal. Concentration} * (\text{Low Next Dil. Ratio} / \text{Cal Nbr 1 Dilution}) = \text{Cal. Concentration} * (10 / 200)$  )
- (\*4) Equal to maximum measuring range multiplied by the re-dilution factor of the primary dilution  
 ( Calculated as:  $\text{Max. Meas. Range} * (\text{Max Next Dil. Ratio} / \text{Sample Dilution}) = \text{Max. Meas. Range} / (120 / 30) =$   
 $= \text{Cal. Concentration} * (\text{Max Next Dil. Ratio} / \text{Cal Nbr 5 Dilution}) = \text{Cal. Concentration} * (120 / 30)$  )
- (\*5) We recommend to disable the automatic control of the calibration interval, and re-calibrate when a new batch of reagents is used, or when the QC established procedures do not give the expected results. Nevertheless, if desired, the user can always define a calibration repeat time in order to be notified of a due calibration when the defined time is elapsed.