3diag - CysC - 800

REF TD-42601

LOT **601-16** 2026-04**(EN)** Cystatin C (serum+urine) - for *IMMAGE*® 800**INSTRUCTIONS FOR USE****(ES)** Cistatina C (suero+orina) - para *IMMAGE*® 800**INSTRUCCIONES DE USO****(PT)** Cistatina C (soro+urina) - para *IMMAGE*® 800**INSTRUÇÕES DE UTILIZAÇÃO**<https://www.3diag.com/006>for IFU, Scan or follow link,
and select LOT*IMMAGE*® is a registered trademark of Beckman Coulter, Inc., Fullerton, CA

TRIMERO Diagnostics, SL

c. València 558, 4t 2a - 08026 Barcelona (Spain)
☎ +34 93 244 86 79 - www.3diag.com

INSTRUCTIONS FOR USE

Reagents for professional use,
for *In Vitro* use only in clinical laboratory (IVD)

3diag - CysC - 800

**Cystatin C (serum+urine)
for IMMAGE® 800**

REF TD-42601

(Product included in **REF** TD-42600)

USO PREVISTO

Quantitative determination of Cystatin C (CysC), in human serum and urine, on **IMMAGE® 800** Immunochemical Systems of **Beckman Coulter** (**IMMAGE®** is a registered trademark of **Beckman Coulter, Inc.**, Fullerton, CA).

PRINCIPLE OF THE METHOD

The specific antibodies (Ab) of the reagent, bound to polystyrene particles, when combined with the antigens (Ag) of the patient sample, form insoluble compounds causing a change in the absorbance and dispersion of the light, proportional to the antigen concentration, which can be quantified by turbidimetric (TIA) or nephelometric (NIA) method, by comparison with calibrators of known concentration.

CONTENTS - COMPOSITION - PREPARATION

- **UDR IMMAGE®** Reagent Cartridge: **REAG UDR CysC**
REF TD-42601-UDR ▼ 100 test containing:
 - Antiserum Reagent: **CONT** 4.5 ml (Compartment A)
Anti-human CysC antibodies bound to polystyrene particles.
 - Accelerator Reagent: **CONT** 1.5 ml (Compartment B)
Solution of Ag/Ab reaction accelerators.

The reagents are ready for use and require no preparation.

Before each use it is convenient that the reagents are homogenized, shaking them gently avoiding the formation of foam or bubbles.

As a preservative, the reagents contain <0.1% (1 g/l) Sodium Azide (NaN_3).

WARNINGS - PRECAUTIONS

- Sodium Azide is toxic. Even if sodium azide is not harmful at the concentration present in the reagents, take the necessary precautions to avoid accidental ingestion or contact with the eyes.
- Sodium Azide can react with lead or copper to give explosive compounds. For disposal it is recommended to rinse with plenty of running water to avoid accumulation in drains.
- Since the absence of infectious agents can not be proven with absolute certainty, components containing materials of human or animal origin must be handled with caution, as potentially infectious, following the recommended safety standards for biological risk.
- Do not mix components belonging to different lot kits.

- Clinical diagnosis should not be based on the results of a single test, but should always integrate all relevant clinical and laboratory data.

STORAGE - SHELF LIFE

- Store refrigerated at +2...+8°C. Do not freeze, as the functionality of the reagents may be altered.
- Properly stored and unopened, the reagents are stable until the expiration date indicated on the label.
- Once opened, the shelf life of the reagents is at least 4 weeks, provided that after each use they are stored immediately in the original containers, tightly capped and refrigerated at +2...+8°C. When using the reagents, it is advisable to use anti-evaporation rubber stoppers (**REF** 447169, from **Beckman Coulter**). The indicated shelf life must be taken as a guideline given that, obviously, it depends on the particular environmental and use conditions, which may differ from those of the stability studies carried out.

MATERIALS NEEDED, NOT SUPPLIED

- **Beckman Coulter's IMMAGE® 800** Immunochemical System, and accessories: dilution segments, cuvettes, etc..
- **BUF2 IMMAGE®**, from **Beckman Coulter** **REF** 447040
- **DIL1 IMMAGE®**, from **Beckman Coulter** **REF** 447640
- **3diag - CysC - Cal Set** **REF** TD-42602
- **3diag - CysC - Control** **REF** TD-42613
- **3diag - CysC - Control (x3)** **REF** TD-42603

SAMPLES

- Fresh Serum.
Samples with presence of fibrin should be centrifuged.
Do not use hemolyzed, lipemic or contaminated samples.
In bibliography⁽¹⁾ a 7-day stability in refrigerated serum is reported.
- Fresh urine.
Prior to the analysis, the samples should be centrifuged until a clear and transparent supernatant is obtained⁽²⁾.
For the determination of specific proteins, centrifugation of urine samples at 3000⁽³⁾-5000⁽⁴⁾ g for 10 minutes is the standard practice in the laboratory.
In bibliography⁽⁵⁾ a 6-day stability in refrigerated urine is reported.

PROCEDURE

Follow the instructions⁽⁶⁾ for use of the analyzer to program and calibrate **UDR** (User Defined Reagent) assays, with the recommended parameters that are detailed below.

Assay Parameters - Serum/Urine

- Units: mg/l
- Protocol: Non-competitive NIPIA
- Sample or Dilution Volume: Serum: 40 µl
Urine: 20 µl
- Reaction Buffer Volume: 195 µl (use **BUF2**)
- Compartment A Volume: 40 µl
- Compartment B Volume: 10 µl
- Gain: 1
- Calibration Dilution: 1:5 (use **DIL1**)
- Sample Dilution: Serum: 1:5 (use **DIL1**)
Urine: 1:1
- Reaction Time: 2.5 minutes

Samples with concentrations higher than the upper limit of the assay range should be analyzed again, diluted manually with **DIL1** to recover a value close to the midpoint of the measurement range.

Calibration Parameters

- Calibrators: Use the **3diag - CysC - Cal Set**.
- Levels: 6
- Replicates: 2
- Model: Logistic 4 Parameters

The assays must be recalibrated, at least when a new batch of reagents is used.

The use of the calibration update option with a single calibrator level is discouraged.

PERFORMANCES OF THE METHOD

Detailed information on the characteristics and performances of the assays is given in the Technical Reports, available on the website (www.3diag.com) or upon request to the Customer Support Service (✉ support@3diag.com - ☎ +34 93 244 86 79).

QUALITY CONTROL

To monitor performances, it is recommended that internal controls be inserted into each analytical series. It is recommended to use the controls of **3diag - CysC - CONTROL**.

Each laboratory should establish its own quality scheme and corrective actions if controls do not meet the assigned tolerances.

The reagents have been subjected to quality control checks and should react as described in these instructions. Therefore, as a general recommendation, in case the controls do not give the expected reaction, as a precaution all reagents should be considered unreliable until their operation has been checked.

TRACEABILITY

Values are referred to the *European Reference Material for the Cystatin C* (cód.: *ERM-DA471/IFCC*) of the *IRMM (Institute for Reference Materials and Measurements)*.

REFERENCE INTERVALS

It is always advisable for each laboratory to establish its own reference values.


The bibliography reports the following reference values:

• Serum⁽¹⁾:

| | | |
|--------|-----------------|------------------|
| Mens: | Up to 22 years: | not established |
| | 23 - 29 years: | 0.60 - 1.03 mg/l |
| | 30 - 39 years: | 0.64 - 1.12 mg/l |
| | 40 - 49 years: | 0.68 - 1.22 mg/l |
| | 50 - 59 years: | 0.72 - 1.32 mg/l |
| | 60 - 69 years: | 0.77 - 1.42 mg/l |
| | 70 - 79 years: | 0.82 - 1.52 mg/l |
| | Over 79 years: | not established |
| Women: | Up to 22 years: | not established |
| | 23 - 29 years: | 0.57 - 0.90 mg/l |
| | 30 - 39 years: | 0.59 - 0.98 mg/l |
| | 40 - 49 years: | 0.62 - 1.07 mg/l |
| | 50 - 59 years: | 0.64 - 1.17 mg/l |
| | 60 - 69 years: | 0.66 - 1.26 mg/l |
| | 70 - 80 years: | 0.68 - 1.36 mg/l |
| | 81 - 86 years: | 0.70 - 1.45 mg/l |
| | Over 86 years: | not established |

The validity of the transfer of these intervals has been verified by analyzing 164 serum samples of presumably healthy adult patients from the Barcelona area, without stratifying them by age, obtaining results (see table and histogram) consistent with those reported in the literature.

| units | mean | SD | range | 95 percentile | 90 percentile |
|-------|-------|-------|--------------|---------------|---------------|
| mg/l | 0.905 | 0.168 | 0.481 - 1.39 | 0.619 - 1.27 | 0.661 - 1.16 |

- Urine⁽⁷⁾: an excretion greather than 0.28 mg/l can be considered significant.

CLINICAL SIGNIFICANCE

Cystatin C, due to its independence from anthropometric variables, constitutes an excellent marker of the glomerular filtration rate (GFR) with which it correlates inversely. The estimation of GFR using Cystatin C, using appropriate equations, is more sensitive, accurate and early than that made using Creatinine, in particular for GFR between 40 and 80 ml/min/1.73m².

Due to its low molecular weight, urine measurement of Cystatin C can be a sensitive marker of tubular dysfunction, due to any cause or pathology.

SYMBOLS

In addition to the harmonized symbols provided on the European Standard EN 980:2008, in the labels and instructions of use has been used the complementary symbology proposed⁽⁸⁾ by the *EDMA (European Diagnostic Manufacturers Association)*, whose meaning is detailed below.

| | |
|-------------|----------------------------|
| REAG | Reagent |
| UDR | UDR (User Defined Reagent) |
| CysC | Cystatin C |
| CONT | Contents |

BIBLIOGRAPHY

- (1) Mayo Medical Laboratories website (www.mayomedicallaboratories.com), date of consultation: 19th September 2017.
- (2) Morales LJ., Ventura S., Solé E et al. - Comité de Comunicación de la Sociedad Española de Medicina de Laboratorio, SEQC^{ML}: "Muestras de Orina de 24 horas y Orina Reciente para la Medición de las Magnitudes Biológicas Más Comunes", ISBN: 978-84-89975-52-1 (2017).
- (3) "Alpha-1-Microglobulin (A1M) - IMMAGE® Immunochemistry Systems Chemistry Information Sheet", © Copyright 2017 Beckman Coulter, Inc..
- (4) Bergón Jiménez E., Bergón Sendín M.: "Uso del cociente cadenas kappa/cadenas lambda en orina para el estudio de la proteína de Bence Jones", Química Clínica 1999; 18 (5) 266-270.
- (5) Pacific Biomarkers website (www.pacbio.com/biomarker/assay-detail/221/), date of consultation: 30th July 2019.
- (6) "IMMAGE® Immunochemical System user's manual". IMMAGE® is a registered trademark of Beckman Coulter, Inc., Fullerton, CA.
- (7) Herget-Rosenthal et al.: "Measurement of urinary Cystatin C by particle-enhanced nephelometric immunoassay: precisión, interferences, stability and reference range", Ann Clin Biochem 2004;41: 111-118.
- (8) EDMA Labelling Task Force: "EDMA Symbols for IVD Reagents and Components - Revision, October 2009".

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
NOTE #1: Applications available for **Alinity c**, **Architect c**, **ADVIA®** series, **AU®** series, **cobas®** series, and other analyzers by request

IMAGE® and AU® are registered trademarks of Beckman Coulter, Inc, Fullerton, CA.
 BN™ and ADVIA® are registered trademarks of Siemens Healthcare Diagnostics Products GmbH, Marburg, Deutschland.
 Optilite® is a registered trademark of The Binding Site Group Ltd, Birmingham, U.K.
 Alinity and Architect related brand marks are registered trademarks of Abbott Laboratories, Abbott Park, Illinois, USA
 cobas® and related brand marks are trademarks of Roche Diagnostics Ltd, Roetkreuz, Switzerland

CONTROLS INCLUDED IN  **Accuracy 365** Online Quality Control Further information and registration at www.accuracy365.com

Accuracy 365 aims to monitor and compare, between laboratories, the results obtained in the Internal Quality Control processes of end-users of specific protein control materials manufactured by **TRIMERO Diagnostics**.

Main features:

- 100 % free and lifetime.
- Cloud application, Available 24 hours, 365 days per year.
- Intuitive interface and very simple to use.
- Calculation and representation of statistical data in real time.
- Anonymous participation, data are only accessible by its owner.