IVD

**3diag - HPX - TIA**

REF TD-42631

LOT **631-02**

⌚ 2023-11

**EN** Hemopexin - for Turbidimetry  
**INSTRUCTIONS FOR USE**

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<https://www.3diag.com/007>

for IFU, Scan or follow link,  
and select LOT

**TRIMERO Diagnostics, SL**

c. València 558, 4t 2a - 08026 Barcelona (Spain)  
☎ +34 93 244 86 79 - [www.3diag.com](http://www.3diag.com)

## INSTRUCTIONS FOR USE

Reagents for professional use,  
for *In Vitro* use only in clinical laboratory (IVD)

## 3diag - HPX - TIA

Hemopexin - for Turbidimetry

**REF TD-42631**

(Product included in Kits, not sold separately)

### INTENDED USE

Quantitative determination of Hemopexin (HPX), in human serum, by turbidimetric method in automatic Clinical Chemistry Analyzers.

### PRINCIPLE OF THE METHOD

The specific antibodies (Ab) of the reagent, when combined with the antigens (Ag) of the patient sample, form insoluble compounds causing a change in the absorbance and dispersion of the light, proportional to the antigen concentration, which can be quantified by turbidimetric (TIA) or nephelometric (NIA) method, by comparison with calibrators of known concentration.

### CONTENTS - COMPOSITION - PREPARATION

- Antiserum Reagent: **REAG Ab HPX**  
**REF** TD-42633-RA    ▽ 100 test<sup>(\*)</sup> - 6 ml  
Solution of anti-human HPX antibodies.
- Reaction Buffer: **BUF HPX**  
**REF** TD-42633-BF    ▽ 100 test<sup>(\*)</sup> - 20 ml  
PBS Buffer, with PEG.

Note (\*1): with the recommended general assay parameters.

The reagents are ready for use and require no preparation.

Before each use it is convenient that the reagents are homogenized, shaking them gently avoiding the formation of foam or bubbles.

As a preservative, the reagents contain <0.1% (1 g/l) Sodium Azide (NaN<sub>3</sub>).

### WARNINGS - PRECAUTIONS

- Sodium Azide is toxic. Even if sodium azide is not harmful at the concentration present in the reagents, take the necessary precautions to avoid accidental ingestion or contact with the eyes.
- Sodium Azide can react with lead or copper to give explosive compounds. For disposal it is recommended to rinse with plenty of running water to avoid accumulation in drains.
- Since the absence of infectious agents can not be proven with absolute certainty, components containing materials of human or animal origin must be handled with caution, as potentially infectious, following the recommended safety standards for biological risk.
- Do not mix components belonging to different lot kits.
- Clinical diagnosis should not be based on the results of a single test, but should always integrate all relevant clinical and laboratory data.

### STORAGE - SHELF LIFE

- Store refrigerated at +2...+8°C. Do not freeze, as the functionality of the reagents may be altered.
- Properly stored and unopened, the reagents are stable until the expiration date indicated on the label.
- Once opened, the shelf life of the reagents is at least 4 weeks, provided that after each use they are stored immediately in the original containers, tightly capped and refrigerated at +2...+8°C. This information should be taken as a guideline given that, obviously, the shelf life depends on the particular environmental and use conditions, which may differ from those of the stability studies carried out.

### MATERIALS NEEDED, NOT SUPPLIED

- Automatic Clinical Chemistry Analyzer, capable of running photometric assays at 340 nm, and accessories: reagent containers, cuvettes, etc..
- The reagents are intended to be used in conjunction with the Calibrators and Controls included in the Kits.

### SAMPLES

Fresh Serum.

Samples with presence of fibrin should be centrifuged.

Do not use hemolyzed, lipemic or contaminated samples.

### PROCEDURE

If necessary, carefully transfer the reagents to the containers used by the analyzer, preventing leakage and foaming or bubbles.

To program and calibrate assays, follow the instructions for use of the analyzer used, with the recommended general parameters that are detailed below. Please, contact the Customer Support Service (✉ [support@3diag.com](mailto:support@3diag.com) - ☎ +34 93 244 86 79) for further information about applications to specific analyzers.

### Assay Parameters

- ① Dispense and mix:
  - Sample/Control: 10 µl (diluted 1:10)
  - Calibrator: 10 µl (diluted 1:10)
  - BUF HPX** 180 µl
- ② Incubate a fixed time between 1 and 5 minutes
- ③ Dispense and mix:
  - REAG Ab HPX** 60 µl
- ④ Read absorbance A1 (Blank) at 340 nm
- ⑤ Incubate a fixed time of about 5 minutes
- ⑥ Read absorbance A2 (End Point) at 340 nm
- ⑦ Interpolate the absorbance increment (A2-A1) of the samples and controls in the curve obtained with the calibrators.
- ⑧ Samples with concentrations higher than the upper limit of the assay range should be analyzed again, diluted manually with Physiological Solution, or by programming a larger sample dilution in the analyzer, to recover a value close to the midpoint of the measurement range.

### Calibration Parameters

- Use the **3diag - HPX - CAL SET**, included in the Kit.
- It is recommended to use Physiological Solution as diluent.
- If the analyzer allows it, it is recommended to program two replicates of each calibration point.
- The calibrations are Non-linear. For the calculation it is recommended to use a 3<sup>rd</sup> Order Polynomial, a Logit or a Polylogal adjustment.

The assay must be recalibrated, at least when a new batch of reagents is used or when its parameterization is changed.

In some analyzers, in order to process the calibrators, it may be necessary to deactivate the clot detection system.

### PERFORMANCES OF THE METHOD

Detailed information on the characteristics and performances of the assay is given in the Technical Report, available on the website ([www.3diag.com](http://www.3diag.com)) or upon request to the Customer Support Service (✉ [support@3diag.com](mailto:support@3diag.com) - ☎ +34 93 244 86 79).

## QUALITY CONTROL

To monitor performances, it is recommended that internal controls be inserted into each analytical series. It is recommended to use the controls of the **3diag - HPX - CONTROL**, included in the Kit.

In some analyzers, in order to process the controls, it may be necessary to deactivate the clot detection system.

Each laboratory should establish its own quality scheme and corrective actions if controls do not meet the assigned tolerances.

The reagents have been subjected to quality control checks and should react as described in these instructions. Therefore, as a general recommendation, in case the controls do not give the expected reaction, as a precaution all reagents should be considered unreliable until their operation has been checked.

## TRACEABILITY

Values in IU/ml are referred to the *1st British Standard for Human Serum Proteins (other than Immunoglobulins) (code: 74/520) (National Institute for Biological Standards and Controls (a World Health Organization (WHO) Laboratory for Biological Standards), NIBSC).*

For the values in mg/dl, given that certified reference materials in these units are not available, the values are referred to internal standards of highly purified proteins. Traceability is ensured by measuring the HPX in the 74/520 standard.

## REFERENCE INTERVALS

It is always advisable for each laboratory to establish its own reference values.

The bibliography<sup>(1)</sup> reports reference values of between 80 to 100 mg/dl, equivalent to about 106 to 133 U/ml.

## CLINICAL SIGNIFICANCE

The serum Hemopexin concentration is a reflection of the amount of Heme released in the blood. After the depletion of Haptoglobin, in severe hemolysis episodes Hemopexin intervenes, which causes its decrease, which can lead to depletion, indicative of the release of large quantities of the Heme group, thus constituting one of the diagnostic indicators of Hemolytic Anemia and a reflection of the severity of hemolysis.

The decrease in Hemopexin is normally associated with the decrease in Haptoglobin except in Thalassemia, in which there is no significant reduction in Haptoglobin.

In pregnancy, Hemopexin levels increase by about 50%.

## SYMBOLS

In addition to the harmonized symbols provided on the European Standard EN 980:2008, in the labels and instructions of use has been used the complementary symbology proposed<sup>(2)</sup> by the *EDMA (European Diagnostic Manufacturers Association)*, whose meaning is detailed below.

**REAG** Reagent

**Ab** Antibody / Antiserum

**BUF** Buffer

**HPX** Hemopexin

## BIBLIOGRAPHY

(1) "Nomenclator de Laboratorio Clínico AEFA/AEBM" - ISBN: 84-4860117-3.

(2) EDMA Labelling Task Force: "EDMA Symbols for IVD Reagents and Components - Revision, October 2009".

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**NOTE #1:** Applications available for **Optilite®**, **Alinity c**, **Architect c**, **ADVIA®** series, **AU®** series, **cobas®** series, and other analyzers by request

**IMAGE®** and **AU®** are registered trademarks of Beckman Coulter, Inc, Fullerton, CA.  
**BN™** and **ADVIA®** are registered trademarks of Siemens Healthcare Diagnostics Products GmbH, Marburg, Deutschland.  
**Optilite®** is a registered trademark of The Binding Site Group Ltd, Birmingham, U.K.  
**Alinity** and **Architect** related brand marks are registered trademarks of Abbott Laboratories, Abbott Park, Illinois, USA  
**cobas®** and related brand marks are trademarks of Roche Diagnostics Ltd, Rottkreuz, Switzerland

**ALL PROTEINS INCLUDED IN**  **Accuracy 365** Online Quality Control Further information and registration at [www.accuracy365.com](http://www.accuracy365.com)

**Accuracy 365** aims to monitor and compare, between laboratories, the results obtained in the Internal Quality Control processes of end-users of specific protein control materials manufactured by **TRIMERO Diagnostics**.

### Main features:

- 100 % free and lifetime.
- Cloud application, Available 24 hours, 365 days per year.
- Intuitive interface and very simple to use.
- Calculation and representation of statistical data in real time.
- Anonymous participation, data are only accessible by its owner.