

3diag - HPX - 800

Hemopexin - for *IMMAGE*[®] 800

REF TD-42621

IVD CE

ANNEX to Instruction for Use:

Application for *IMMAGE*[®] and *IMMAGE*[®] 800 Immunochemical Systems

<u>PROTOCOL INFORMATION</u>			
Reagent Lot Number:	(^{*1})	Units:	mg/dL or U/mL
Cartridge Lot Number:	(^{*1})	Protocol:	Non-competitive Nephelometry
Cartridge Serial Number:	(^{*1})	AGXS Limit:	not defined
Reagent Expiry date:	(^{*1})		
Sample or Dilution Volume:	24 µL	Diluent:	DIL1 (^{*2})
Reaction Buffer Volume:	195 µL	Buffer:	BUF1 (^{*3})
Calibration Dilution:	1 : 24	Gain:	3
Compartment A Volume:	35 µL	Sample Dilution:	1 : 24
Compartment B Volume:	5 µL	Reaction Time:	2.0 minutes
<u>CALIBRATION INFORMATION</u>			
Levels:	6	Replicates:	2
Update Level:	4 (^{*4})	Replicates:	2 (^{*4})
Model:	Logistic 4 parameters or 3rd order Polynomial (select the best adjustment)		
<u>NOTES</u>			
<p>(^{*1}): See the labels of the Reagent Cartridge. (^{*2}): <i>Beckman Coulter</i> P/N no. 447640. For the <i>IMMAGE</i>[®] it has to be defined as DIL10. (^{*3}): <i>Beckman Coulter</i> P/N no. 447650. For the <i>IMMAGE</i>[®] it has to be defined as BUF10. (^{*4}): The use of the calibration update, with only one calibrator level, is discouraged.</p>			