**3diag - RBP - TIA**

REF TD-42671

LOT **671-23**

 2025-02

UDI

(01) 08434477309052  
(10) 671-23  
(17) 250228  
(240) TD-42671

- (DE)** Retinol-bindendes Protein (Serum+Urin) - für Turbidimetrie  
**GEBRAUCHSANWEISUNGEN**
- (EN)** Retinol Binding Protein (serum+urine) - for Turbidimetry  
**INSTRUCTIONS FOR USE**
- (ES)** Proteína Fijadora del Retinol (suero+orina) - para Turbidimetría  
**INSTRUCCIONES DE USO**
- (FR)** Protéine de Liaison au Rétinol (sérum+urine) - pour Turbidimétrie  
**INSTRUCTIONS D'UTILISATION**
- (IT)** Proteina Legante il Retinolo (siero+urina) - per Turbidimetria  
**ISTRUZIONI PER L'USO**

<https://www.3diag.com/009>

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#### BIBLIOGRAPHIE

- (1) Webseite von *Mayo Medical Laboratories* ([www.mayomedicallaboratories.com](http://www.mayomedicallaboratories.com)), Abrufdatum: 18. April 2017.
- (2) Miura N. et al.: "Evaluation of the reference range of retinol-binding protein (RBP) levels by the latex turbidimetric immunoassay", *Rinsho Byori* 2009; 57(3): 195-9.
- (3) Siemens Healthcare Diagnostics Products GmbH: "Sistema BN<sup>®</sup> II - Protocolos de Ensayo - Versión 2.4", 2009/03.
- (4) Webseite von *ARUP Laboratories* ([www.aruplab.com](http://www.aruplab.com)), Abrufdatum: 24. Januar 2019.
- (5) Forga L. et al.: "Low Serum Levels of Prealbumin, Retinol Binding Protein, and Retinol Are Frequent in Adult Type 1 Diabetic Patients", *Journal of Diabetes Research* 2016; Article ID 2532108.
- (6) Bernard AM. et al.: "Assessment of Urinary Retinol-Binding Protein as an Index of Proximal Tubular Injury", *Clin. Chem.* 33/6, 775-779 (1987).
- (7) EDMA Labelling Task Force: "EDMA Symbols for IVD Reagents and Components - Revision, October 2009".

#### REVISIONSDATUM DES TEXTES

25. Juli 2022.

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## INSTRUCTIONS FOR USE

Reagents for professional use,  
for *In Vitro* use only in clinical laboratory (IVD)

# 3diag - RBP - TIA

Retinol Binding Protein  
for Turbidimetry

**REF TD-42671**

### INTENDED USE

Quantitative determination of Retinol Binding Protein (RBP), in human serum and urine, by turbidimetric method in automatic Clinical Chemistry Analyzers.

### PRINCIPLE OF THE METHOD

The specific antibodies (Ab) of the reagent, bound to polystyrene particles, when combined with the antigens (Ag) of the patient sample, form insoluble compounds causing a change in the absorbance and dispersion of the light, proportional to the antigen concentration, which can be quantified by turbidimetric (TIA) or nephelometric (NIA) method, by comparison with calibrators of known concentration.

### CONTENTS - COMPOSITION - PREPARATION

- Antiserum Reagent: **REAG Ab RBP**  
**REF** TD-42671-RA ▽ 100 test <sup>(\*)</sup> - 3.5 ml  
Anti-RBP antibodies, bound to polystyrene particles.
- Reaction Buffer: **BUF RBP**  
**REF** TD-42671-BF ▽ 100 test <sup>(\*)</sup> - 20 ml  
TRIS Buffer, with PEG.

Note (\*1): with the recommended general assay parameters.

As a preservative, the reagents contain <0.1% (1 g/l) Sodium Azide (NaN<sub>3</sub>).

The reagents are ready for use and require no preparation.

Before each use it is convenient that the reagents are homogenized, shaking them gently avoiding the formation of foam or bubbles.

### WARNINGS - PRECAUTIONS

- Sodium Azide is toxic. Even if sodium azide is not harmful at the concentration present in the reagents, take the necessary precautions to avoid accidental ingestion or contact with the eyes.
- Sodium Azide can react with lead or copper to give explosive compounds. For disposal it is recommended to rinse with plenty of running water to avoid accumulation in drains.
- Since the absence of infectious agents can not be proven with absolute certainty, components containing materials of human or animal origin must be handled with caution, as potentially infectious, following the recommended safety standards for biological risk.
- Do not mix components belonging to different lot kits.
- Clinical diagnosis should not be based on the results of a single test, but should always integrate all relevant clinical and laboratory data.

### STORAGE - SHELF LIFE

- Store refrigerated at +2...+8°C. Do not freeze, as the functionality of the reagents may be altered.
- Properly stored and unopened, the reagents are stable until the expiration date indicated on the label.
- Once opened, the shelf life of the reagents is at least 4 weeks, provided that after each use they are stored immediately in the original containers, tightly capped and refrigerated at +2...+8°C. This information should be taken as a guideline given that, obviously, the shelf life depends on the particular environmental and use conditions, which may differ from those of the stability studies carried out.

### MATERIALS NEEDED, NOT SUPPLIED

- Automatic Clinical Chemistry Analyzer, capable of running photometric assays at 540...600 nm, and accessories: reagent containers, cuvettes, etc..

- 3diag - S-RBP - CAL SET** **REF** TD-42662-S
- 3diag - RBP - CAL** **REF** TD-42672
- 3diag - RBP - CONTROL** **REF** TD-42663
- 3diag - S-RBP - CONTROL** **REF** TD-42673-S

### SAMPLES

- Fresh Serum.  
Samples with presence of fibrin should be centrifuged.  
Do not use hemolyzed, lipemic or contaminated samples.

- Fresh Urine.

Samples should be centrifuged.

In bibliography<sup>(1)</sup> it is reported the following stability :

- Serum: Refrigerated: 7 days - Frozen: 90 days
- Urine: Refrigerated/Frozen: 7 days

### PROCEDURE

If necessary, carefully transfer the reagents to the containers used by the analyzer, preventing leakage and foaming or bubbles. To program and calibrate assays, follow the instructions for use of the analyzer used, with the recommended general parameters that are detailed below. Please, contact the Customer Support Service (✉ [support@3diag.com](mailto:support@3diag.com) - ☎ +34 93 244 86 79) for further information about applications to specific analyzers.

### Assay Parameters - Serum/Urine

- ①Dispense and mix:
  - \* **Serum:** Sample/Control: 15 ul (diluted 1:30)  
Calibrator: 15 ul (see dilution factor in Calibration Parameters)
  - \* **Urine:** Sample/Control: 15 ul (diluted 1:5)  
Calibrator: 15 ul (see dilution factor in Calibration Parameters)
  - \* **BUF RBP** 180 ul
- ②Incubate a fixed time between 1 and 5 minutes
- ③Dispense and mix:
  - \* **REAG Ab RBP** 30 ul
- ④Read absorbance A1 (Blank) at 540...600 nm
- ⑤Incubate a fixed time of about 5 minutes
- ⑥Read absorbance A2 (Final Point) at 540...600 nm
- ⑦Interpolate the absorbance increment (A2-A1) of the samples and controls in the curve obtained with the calibrators
- ⑧Samples with concentrations higher than the upper limit of the assay range should be analyzed again, diluted manually with saline, or by programming a larger sample dilution in the analyzer, to recover a value close to the midpoint of the measurement range

### Calibration Parameters - Serum/Urine

- Serum:**  
If using the **3diag - S-RBP - CAL SET** program a dilution factor of 1:5.  
If using the **3diag - RBP - CAL** program the following dilutions: 1:5, 1:7, 1:15, 1:30, 1:60 and 1:120.
- Urine:**  
If using the **3diag - S-RBP - CAL SET** program a dilution factor of 1:5.  
If using the **3diag - RBP - CAL** program the following dilutions: 1:5, 1:10, 1:20, 1:40, 1:120 and 1:700 (alternatively, saline can be used as the lowest calibration point - zero value).
- It is recommended to use saline as diluent.

- If the analyzer allows it, it is recommended to program two replicates of each calibration point.
- The calibrations are Non-linear. For the calculation it is recommended to use a 3<sup>rd</sup> Order Polynomial, a Logit or a Polygonal adjustment.
- The assay must be recalibrated, at least when a new batch of reagents is used or when its parameterization is changed.

#### PERFORMANCES OF THE METHOD

Detailed information on the characteristics and performances of the assays is given in the Technical Reports, available on the website ([www.3diag.com](http://www.3diag.com)) or upon request to the Customer Support Service (✉ [support@3diag.com](mailto:support@3diag.com) - ☎ +34 93 244 86 79).

#### QUALITY CONTROL

To monitor performances, it is recommended that internal controls be inserted into each analytical series. It is recommended to use the controls of the **3diag - RBP - CONTROL** or the **3diag - S-RBP - CONTROL**.

Each laboratory should establish its own quality scheme and corrective actions if controls do not meet the assigned tolerances.

The reagents have been subjected to quality control checks and should react as described in these instructions. Therefore, as a general recommendation, in case the controls do not give the expected reaction, as a precaution all reagents should be considered unreliable until their operation has been checked.

#### TRACEABILITY

Given that certified reference materials are not available, values are referred to internal standards based on highly purified proteins.

Values are also referred to commercial calibrators *N Protein Standard SL* (REF: *OQIM13*, LOT: *083643*) from *Siemens Healthcare*.

Traceability is ensured by measuring the RBP in the *European Reference Material ERM-DA470k/IFCC* (*Institute for Reference Materials and Measurements, IRMM*).

#### REFERENCE INTERVALS

It is always advisable for each laboratory to establish its own reference values.


##### Serum

The bibliography reports variable reference values between publications, depending on the method and standardization used and the population analyzed:

- between 1.5 to 6.7 mg / dl <sup>(1)</sup>,
- 2.65 to 6.0 mg/dl for men and 1.9 to 4.61 mg/dl for women <sup>(2)</sup>,
- 3.0 to 6.0 mg/dl <sup>(3)(4)(5)</sup>,
- and 2.0 to 5.0 mg/dl <sup>(6)</sup>.

Analyzing serum samples from 110 presumably healthy adult patients from the Barcelona area, without stratifying them by age or sex, the following results have been obtained (see table and histogram - values based on *Trimerio Diagnostics'* internal standardization):

units	mean	SD	range	95 percentile	90 percentile
mg/dl	3.085	0.624	1.91 - 5.63	2.08 - 4.25	2.15 - 4.18

#### Urine

The bibliography<sup>(1)</sup> reports that an excretion greater than 0.163 mg/24h or 0.130 mg/g-creatinine can be considered significant.

#### CLINICAL SIGNIFICANCE

Serum RBP is elevated in Diabetes Mellitus tipo 2, obesity, metabolic syndrome, cardiovascular disease and chronic kidney disease. Decreased values due to proteinuria, hepatocellular disease, hyperparathyroidism, cystic fibrosis or other causes, can pose a serious threat to transportation of Vitamin A to tissues.

Urinary RBP results a good marker of tubular dysfunction, due to any cause or pathology, given its good sensitivity, stability and high dynamic range (complete failure of tubular reabsorption increases its excretion in a factor of 10<sup>4</sup>-10<sup>5</sup> times the normal values).

#### SYMBOLS

In addition to the harmonized symbols provided on the European Standard EN 980:2008, in the labels and instructions of use has been used the complementary symbology proposed<sup>(7)</sup> by the *EDMA* (*European Diagnostic Manufacturers Association*), whose meaning is detailed below.

<b>REAG</b>	Reagent
<b>Ab</b>	Antibody / Antiserum
<b>BUF</b>	Buffer
<b>RBP</b>	Retinol Binding Protein

#### BIBLIOGRAPHY

- (1) *Mayo Medical Laboratories* website ([www.mayomedicallaboratories.com](http://www.mayomedicallaboratories.com)), date of consultation: 18<sup>th</sup> April 2017.
- (2) Miura N. et al.: "Evaluation of the reference range of retinol-binding protein (RBP) levels by the latex turbidimetric immunoassay", *Rinsho Byori* 2009; 57(3): 195-9.
- (3) *Siemens Healthcare Diagnostics Products GmbH*: "Sistema BN<sup>o</sup> II - Protocolos de Ensayo - Versión 2.4", 2009/03.
- (4) *ARUP Laboratories* website ([www.aruplab.com](http://www.aruplab.com)), date of consultation: 24<sup>th</sup> January 2019.
- (5) Forga L. et al.: "Low Serum Levels of Prealbumin, Retinol Binding Protein, and Retinol Are Frequent in Adult Type 1 Diabetic Patients", *Journal of Diabetes Research* 2016; Article ID 2532108.
- (6) Bernard AM. et al.: "Assessment of Urinary Retinol-Binding Protein as an Index of Proximal Tubular Injury", *Clin. Chem.* 33/6, 775-779 (1987).
- (7) *EDMA Labelling Task Force*: "EDMA Symbols for IVD Reagents and Components - Revision, October 2009".

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#### BIBLIOGRAPHIE

- (1) Page Web de *Mayo Medical Laboratories* ([www.mayomedicallaboratories.com](http://www.mayomedicallaboratories.com)), date consultation: 18/04/2017.
- (2) Miura N. et al.: "Evaluation of the reference range of retinol-binding protein (RBP) levels by the latex turbidimetric immunoassay", *Rinsho Byori* 2009; 57(3): 195-9.
- (3) Siemens Healthcare Diagnostics Products GmbH: "Sistema BN° II - Protocolos de Ensayo - Versión 2.4", 2009/03.
- (4) Page Web de *ARUP Laboratories* ([www.aruplab.com](http://www.aruplab.com)), date consultation: 24/01/2019.
- (5) Forga L. et al.: "Low Serum Levels of Prealbumin, Retinol Binding Protein, and Retinol Are Frequent in Adult Type 1 Diabetic Patients", *Journal of Diabetes Research* 2016; Article ID 2532108.
- (6) Bernard AM. et al.: "Assessment of Urinary Retinol-Binding Protein as an Index of Proximal Tubular Injury", *Clin. Chem.* 33/6, 775-779 (1987).
- (7) EDMA Labelling Task Force: "EDMA Symbols for IVD Reagents and Components - Revision, October 2009".

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## BIBLIOGRAFIA

- (1) Sito Web di **Mayo Medical Laboratories** ([www.mayomedicallaboratories.com](http://www.mayomedicallaboratories.com)), data consultazione: 18 aprile 2017.
- (2) Miura N. et al.: "Evaluation of the reference range of retinol-binding protein (RBP) levels by the latex turbidimetric immunoassay", *Rinsho Byori* 2009; 57(3): 195-9.
- (3) Siemens Healthcare Diagnostics Products GmbH: "Sistema BN® II - Protocolos de Ensayo - Versión 2.4", 2009/03.
- (4) Sito Web di **ARUP Laboratories** ([www.aruplab.com](http://www.aruplab.com)), data consultazione: 24 gennaio 2019.
- (5) Forga L. et al.: "Low Serum Levels of Prealbumin, Retinol Binding Protein, and Retinol Are Frequent in Adult Type 1 Diabetic Patients", *Journal of Diabetes Research* 2016; Article ID 2532108.
- (6) Bernard AM. et al.: "Assessment of Urinary Retinol-Binding Protein as an Index of Proximal Tubular Injury", *Clin. Chem.* 33/6, 775-779 (1987).
- (7) EDMA Labelling Task Force: "EDMA Symbols for IVD Reagents and Components - Revision, October 2009".

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25 luglio 2022.

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unità	media	DS	intervallo	percentile 95	percentile 90
mg/dl	3.085	0.624	1.91 - 5.63	2.08 - 4.25	2.15 - 4.18

## CATALOGUE of Available Proteins - Turbidimetry (TIA) / Nephelometry (NIA)

Protein Description	Traced to	Sample	Beckman C. IMAGE®	Siemens H. BN™ series	Turbidimetry see NOTE #1
 Free Light Chains - Serum	ERM-DA470k	Serum	Available	-	Available
 Free Light Chains - Urine	ERM-DA470k	Urine	Available	Available	Available
Kappa (B&F) Light Chains	ERM-DA470k	Serum & Urine	-	-	Available
Lambda (B&F) Light Chains	ERM-DA470k	Serum & Urine	-	-	Available
Beta-2 Microglobulin	WHO B2M	Serum & Urine	Available	-	Available
C1q Complement	WHO W1032	Serum	Available	Available	Available
C5 Complement	WHO W1032	Serum	Available	Available	Available
C1 (Esterase) Inhibitor	Internal Standard	Serum	Available	-	Available
Factor B (C3 Proactivator)	WHO W1032	Serum	Available	Available	Available
Cystatin C	ERM-DA471	Serum & Urine	Available	only in Urine	Available
Hemopexin	NIBSC 74/520	Serum	Available	-	Available
IgD Immunoglobulins	NIBSC 67/037	Serum	Available	Available	Available
Retinol Binding Protein	Internal Standard	Serum & Urine	Available	only in Urine	Available
Soluble Transferrin Receptor	WHO 07/202	Serum	Available	-	Available
IgE Immunoglobulins	WHO 75/502	Serum	-	-	Available
Alpha-1 Microglobulin	Internal Standard	Urine	-	-	Available
Serum Amyloid A	WHO 92/680	Serum	Available	-	Available

**NOTE #1:** Applications available for **Optilite®**, **Alinity c**, **Architect c**, **ADVIA®** series, **AU®** series, **cobas®** series, and other analyzers by request

**IMAGE®** and **AU®** are registered trademarks of Beckman Coulter, Inc, Fullerton, CA.  
**BN™** and **ADVIA®** are registered trademarks of Siemens Healthcare Diagnostics Products GmbH, Marburg, Deutschland.  
**Optilite®** is a registered trademark of The Binding Site Group Ltd, Birmingham, U.K.  
**Alinity** and **Architect** related brand marks are registered trademarks of Abbott Laboratories, Abbott Park, Illinois, USA  
**cobas®** and related brand marks are trademarks of Roche Diagnostics Ltd, Rottkreuz, Switzerland

**ALL PROTEINS INCLUDED IN**  **Accuracy 365** Online Quality Control Further information and registration at [www.accuracy365.com](http://www.accuracy365.com)

**Accuracy 365** aims to monitor and compare, between laboratories, the results obtained in the Internal Quality Control processes of end-users of specific protein control materials manufactured by **TRIMERO Diagnostics**.

### Main features:

- 100 % free and lifetime.
- Cloud application, Available 24 hours, 365 days per year.
- Intuitive interface and very simple to use.
- Calculation and representation of statistical data in real time.
- Anonymous participation, data are only accessible by its owner.