**3diag - SAA - 800**

REF TD-42881

LOT **881-36**

2026-03

EN Serum Amyloid A - for *IMMAGE*<sup>®</sup> 800  
**INSTRUCTIONS FOR USE**

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<https://www.3diag.com/018>  
for IFU, Scan or follow link,  
and select LOT

*IMMAGE*<sup>®</sup> is a registered trademark of Beckman Coulter, Inc., Fullerton, CA

**TRIMERO Diagnostics, SL**

c. València 558, 4t 2a - 08026 Barcelona (Spain)  
 ☎ +34 93 244 86 79 - [www.3diag.com](http://www.3diag.com)

## INSTRUCTIONS FOR USE

Reagents for professional use,  
 for *In Vitro* use only in clinical laboratory (IVD)

# 3diag - SAA - 800

**Serum Amyloid A  
 for IMMAGE® 800**

**REF TD-42881**

(Product included in REF TD-42880, not sold separately)

### INTENDED USE

Quantitative determination of Serum Amyloid A (SAA), in human serum, on IMMAGE® 800 Immunochemical Systems of Beckman Coulter (IMMAGE® is a registered trademark of Beckman Coulter, Inc., Fullerton, CA).

### PRINCIPLE OF THE METHOD

The specific antibodies (Ab) of the reagent, bound to polystyrene particles, when combined with the antigens (Ag) of the patient sample, form insoluble compounds causing a change in the absorbance and dispersion of the light, proportional to the antigen concentration, which can be quantified by turbidimetric (TIA) or nephelometric (NIA) method, by comparison with calibrators of known concentration.

### CONTENTS - COMPOSITION - PREPARATION

- UDR IMMAGE® Reagent Cartridge: **REAG UDR SAA**  
 REF TD-42881-UDR ▼ 100 test containing:
  - Antiserum Reagent: **CONT** 2.5 ml (Compartment A)  
 Solution of anti-SAA antibodies, bound to polystyrene particles.
  - Accelerator Reagent: **CONT** 2.5 ml (Compartment B)  
 Solution of Ag/Ab reaction accelerators.
- Reaction Buffer: **BUF SAA**  
 REF TD-42881-BUF **CONT** 32 ml  
 TRIS Buffer, with PEG.

The reagents are ready for use and require no preparation. Before each use it is convenient that the reagents are homogenized, shaking them gently avoiding the formation of foam or bubbles.

As a preservative, the reagents contain <0.1% (1 g/l) Sodium Azide (NaN<sub>3</sub>).

### WARNINGS - PRECAUTIONS

- Sodium Azide is toxic. Even if sodium azide is not harmful at the concentration present in the reagents, take the necessary precautions to avoid accidental ingestion or contact with the eyes.
- Sodium Azide can react with lead or copper to give explosive compounds. For disposal it is recommended to rinse with plenty of running water to avoid accumulation in drains.

- Since the absence of infectious agents can not be proven with absolute certainty, components containing materials of human or animal origin must be handled with caution, as potentially infectious, following the recommended safety standards for biological risk.
- Do not mix components belonging to different lot kits.
- Clinical diagnosis should not be based on the results of a single test, but should always integrate all relevant clinical and laboratory data.

### STORAGE - SHELF LIFE

- Store refrigerated at +2...+8°C. Do not freeze, as the functionality of the reagents may be altered.
- Properly stored and unopened, the reagents are stable until the expiration date indicated on the label.
- Once opened, the shelf life of the reagents is at least 4 weeks, provided that after each use they are stored immediately in the original containers, tightly capped and refrigerated at +2...+8°C. When using the reagents, it is advisable to use anti-evaporation rubber stoppers (REF 447169, from Beckman Coulter). The indicated shelf life must be taken as a guideline given that, obviously, it depends on the particular environmental and use conditions, which may differ from those of the stability studies carried out.

### MATERIALS NEEDED, NOT SUPPLIED

- Beckman Coulter's IMMAGE® 800 Immunochemical System, and accessories: dilution segments, cuvettes, etc..
- DIL1 IMMAGE®, from Beckman Coulter REF 447640
- 3diag - SAA - Cal Set REF TD-42882
- 3diag - SAA - Control REF TD-42883

### SAMPLES

Fresh Serum.

#### Bring samples to room temperature before analysis.

Samples with presence of fibrin should be centrifuged. Do not use hemolyzed, lipemic or contaminated samples.

In bibliography<sup>(1)</sup> it is reported the following stability in serum:

- Refrigerated: 8 days
- Freezed (-20°C): 1 year
- Freezed (-70°C): 4.5 years

### PROCEDURE

Follow the instructions<sup>(2)</sup> for use of the analyzer to program and calibrate UDR (User Defined Reagent) assays, with the recommended parameters that are detailed below.

#### Assay Parameters

- Units: mg/l
- Protocol: Non-competitive Nephelometry
- Sample or Dilution Volume: 5 µl
- Reaction Buffer Volume: 200 µl  
 (use **BUF SAA**, defined as BUF10 to BUF15)
- Compartment A Volume: 16 µl
- Compartment B Volume: 20 µl
- Gain: 1
- Calibration Dilution: 1:1 (neat)
- Sample Dilution: 1:1 (neat)
- Reaction Time: 5.0 minutes

Samples with concentrations higher than the upper limit of the assay range should be analyzed again, diluted manually with DIL1, or by programming a larger sample dilution in the analyzer, to recover a value close to the midpoint of the measurement range.

#### Calibration Parameters

- Calibrators: Use the 3diag - SAA - Cal Set
  - Levels: 6 - Replicates: 2 - Model: Logistic 4 Parameters
- The assay must be recalibrated, at least when a new batch of reagents is used, or when the internal quality control established procedures do not give the expected results. The use of the calibration update option with a single calibrator level is discouraged.

## PERFORMANCES OF THE METHOD

Detailed information on the characteristics and performances of the assay is given in the Technical Report, available on the website ([www.3diag.com](http://www.3diag.com)) or upon request to the Customer Support Service (✉ [support@3diag.com](mailto:support@3diag.com) - ☎ +34 93 244 86 79).

## QUALITY CONTROL

To monitor performances, it is recommended that internal controls be inserted into each analytical series. It is recommended to use the controls of **3diag - SAA - CONTROL**.

Each laboratory should establish its own quality scheme and corrective actions if controls do not meet the assigned tolerances.

The reagents have been subjected to quality control checks and should react as described in these instructions. Therefore, as a general recommendation, in case the controls do not give the expected reaction, as a precaution all reagents should be considered unreliable until their operation has been checked.

## TRACEABILITY

Values are referred to the *Serum Amyloid A (SAA) 1st International Standard* (NIBSC code: 92/680) of the *WHO (World Health Organisation)*.

## REFERENCE INTERVALS


It is always advisable for each laboratory to establish its own reference values.

In the literature it is reported that, under normal conditions, the SAA circulates at the level of traces (up to 3-10 mg/l, depending on the method used and the population analyzed); for example:

- up to 5 mg/l<sup>(1)</sup> in ELISA,
- up to 10 mg/l<sup>(3)</sup> (method not specified),
- up to 7.4 mg/l<sup>(4)</sup>, 6.4 mg/l<sup>(5)</sup> and 8.0 mg/l<sup>(6)</sup> in nephelometry.

Analyzing, with turbidimetric method, serum samples from 120 presumably healthy patients from the Barcelona area (discarding those with C-reactive Protein out of the reference interval), the following results (see table and histogram) have been obtained (units are mg/l):


mean	SD	range	95 percentile	90 percentile
3.81	2.20	0.84 - 14.2	up to <b>8.25</b>	up to <b>6.62</b>

In view of the results, a concentration higher than about 8 mg/l can be considered significant.

The validity of the transference of the above intervals to the Rate-Nephelometric method has been verified by a correlation study with 40 samples, obtaining the following results (see table and regression plot):

Correlation Coefficient	Linear Regression	R <sup>2</sup>	N
0.9909	y = 0.9284 x + 1.5405	0.9818	40

## CLINICAL SIGNIFICANCE

Serum Amyloid A (SAA) is an acute-phase protein. During acute events, the rise in SAA levels is one of the most rapid and intense increases of all acute-phase proteins. Only a few hours after inflammatory stimulus, SAA levels can increase by as much as 1000 fold thus making SAA a sensitive marker of inflammatory response. The acute-phase response usually lasts for several days and then the concentration of SAA gradually decreases in the absence of a new stimulus.

Measuring SAA levels may be a useful indicator of the response to therapy and degree of acute and chronic inflammation, due to any inflammatory disorder, such as rheumatoid arthritis, juvenile arthritis, ankylosing spondylitis, familial Mediterranean fever, progressive sclerosis as well as bacterial infections. Secondary amyloidosis may develop as a result of prolonged or repeated inflammatory conditions in which SAA levels remain elevated.

## SYMBOLS

In addition to the harmonized symbols provided on the European Standard EN 980:2008, in the labels and instructions of use has been used the complementary symbology proposed<sup>(7)</sup> by the *EDMA (European Diagnostic Manufacturers Association)*, whose meaning is detailed below.

<b>REAG</b>	Reagent
<b>UDR</b>	UDR (User Defined Reagent)
<b>SAA</b>	Serum Amyloid A
<b>BUF</b>	Buffer
<b>CONT</b>	Contents

## BIBLIOGRAPHY

- (1) *Pacific Biomarkers* website ([www.pacbio.com/biomarker/assay-detail/63/](http://www.pacbio.com/biomarker/assay-detail/63/)), date of consultation: 17<sup>th</sup> November 2017.
- (2) "IMMAGE® Immunochemical System user's manual". IMMAGE® is a registered trademark of Beckman Coulter, Inc., Fullerton, CA.
- (3) B. Targonska-Stepniak and M. Majdan: "Serum Amyloid A as a Marker of Persistent Inflammation and an Indicator of Cardiovascular and Renal Involvement in Patients with Rheumatoid Arthritis" - *Mediators of Inflammation* (Hindawi P.C.) 2014; Art. ID 793628 (<http://dx.doi.org/10.1155/2014/793628>).
- (4) M.L. Seco, L. Borque et al.: "Evaluación de un nuevo método inmunofelométrico para la determinación de proteína amiloide sérica A" - *Química Clínica* 2002; 21 (1) 5-9.
- (5) Siemens Healthcare Diagnostics Products GmbH: "Sistema BN° II - Protocolos de Ensayo - Versión 2.4", 2009/03.
- (6) JY. Wang et al: "Elevated levels of serum amyloid A indicate poor prognosis in patients with esophageal squamous cell carcinoma" - *BMC Cancer* 2012; 12:365 ([www.biomedcentral.com/1471-2407/12/365](http://www.biomedcentral.com/1471-2407/12/365)).
- (7) EDMA Labelling Task Force: "EDMA Symbols for IVD Reagents and Components - Revision, October 2009".

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**NOTE #1:** Applications available for **Alinity c**, **Architect c**, **ADVIA®** series, **AU®** series, **cobas®** series, and other analyzers by request

**IMMAGE®** and **AU®** are registered trademarks of Beckman Coulter, Inc, Fullerton, CA.  
**BN™** and **ADVIA®** are registered trademarks of Siemens Healthcare Diagnostics Products GmbH, Marburg, Deutschland.  
**Optilite®** is a registered trademark of The Binding Site Group Ltd, Birmingham, U.K.  
**Alinity** and **Architect** related brand marks are registered trademarks of Abbott Laboratories, Abbott Park, Illinois, USA  
**cobas®** and related brand marks are trademarks of Roche Diagnostics Ltd, Roetkreuz, Switzerland

**CONTROLS INCLUDED IN**  **Accuracy 365** Online Quality Control Further information and registration at [www.accuracy365.com](http://www.accuracy365.com)

**Accuracy 365** aims to monitor and compare, between laboratories, the results obtained in the Internal Quality Control processes of end-users of specific protein control materials manufactured by **TRIMERO Diagnostics**.

### Main features:

- 100 % free and lifetime.
- Cloud application, Available 24 hours, 365 days per year.
- Intuitive interface and very simple to use.
- Calculation and representation of statistical data in real time.
- Anonymous participation, data are only accessible by its owner.

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- (5) Siemens Healthcare Diagnostics Products GmbH: "Sistema BN® II - Protocolos de Ensayo - Versión 2.4", 2009/03.