3diag - sTfR - ADV(5)

REF TD-42696

LOT **696-14**

⏳ 2024-10

UDI

(01) 08434477310089
(10) 696-14
(17) 241030
(240) 11415057

EN Soluble Transferrin Receptor - for **ADVIA**[®] (500 test)
INSTRUCTIONS FOR USE

ES Receptor Soluble de la Transferrina - para **ADVIA**[®] (500 test)
INSTRUCCIONES DE USO

<https://www.3diag.com/010>

for IFU, Scan or follow link,
and select LOT

ADVIA[®] and its associated brands are registered trademarks of *Siemens Healthcare*

TRIMERO Diagnostics, SL

c. València 558, 4t 2a - 08026 Barcelona (Spain)
☎ +34 93 244 86 79 - www.3diag.com

INSTRUCTIONS FOR USE

Reagents for professional use,
for *In Vitro* use only in clinical laboratory (IVD)

3diag - sTfR - ADV (5)

**Soluble Transferrin Receptor
for ADVIA® (500 test)**

REF TD-42696

INTENDED USE

Quantitative determination of Soluble Transferrin Receptor (sTfR), in human serum, in the ADVIA® Chemistry Systems series clinical chemistry analyzers of Siemens Healthcare (ADVIA® and its associated brands are registered trademarks of Siemens Healthcare).

PRINCIPLE OF THE METHOD

The specific antibodies (Ab) of the reagent, bound to polystyrene particles, when combined with the antigens (Ag) of the patient sample, form insoluble compounds causing a change in the absorbance and dispersion of the light, proportional to the antigen concentration, which can be quantified by turbidimetric (TIA) or nephelometric (NIA) method, by comparison with calibrators of known concentration.

CONTENTS - COMPOSITION - PREPARATION

- Antiserum Reagent: **REAG Ab sTfR**
REF TD-42696-RA ▽ 500 test - 15 ml
Anti-human sTfR antibodies, bound to polystyrene particles.
- Reaction Buffer: **BUF sTfR**
REF TD-42696-BF ▽ 500 test - 55 ml
TRIS buffer, with PEG.

As a preservative, the reagents contain <0.1% (1 g/l) Sodium Azide (NaN₃).

The reagents are ready for use and require no preparation. Before each use it is convenient that the reagents are homogenized, shaking them gently avoiding the formation of foam or bubbles.

WARNINGS - PRECAUTIONS

- Sodium Azide is toxic. Even if sodium azide is not harmful at the concentration present in the reagents, take the necessary precautions to avoid accidental ingestion or contact with the eyes.
- Sodium Azide can react with lead or copper to give explosive compounds. For disposal it is recommended to rinse with plenty of running water to avoid accumulation in drains.
- Since the absence of infectious agents can not be proven with absolute certainty, components containing materials of human or animal origin must be handled with caution, as potentially infectious, following the recommended safety standards for biological risk.
- Do not mix components belonging to different lot kits.
- Clinical diagnosis should not be based on the results of a single test, but should always integrate all relevant clinical and laboratory data.

STORAGE - SHELF LIFE

- Store refrigerated at +2...+8°C. Do not freeze, as the functionality of the reagents may be altered.
- Properly stored and unopened, the reagents are stable until the expiration date indicated on the label.
- Once opened, the shelf life of the reagents is at least 4 weeks, provided that after each use they are stored immediately in the original containers, tightly capped and refrigerated at +2...+8°C. This information should be taken as a guideline given that, obviously, the shelf life depends on the particular environmental and use conditions, which may differ from those of the stability studies carried out.

MATERIALS NEEDED, NOT SUPPLIED

- ADVIA® Chemistry Systems series analyzer of Siemens Healthcare, and accessories: adapters, cuvettes, etc..
- 3diag - sTfR - CAL SET** **REF** TD-42682
- 3diag - sTfR - CONTROL** **REF** TD-42683

SAMPLES

- Fresh Serum.
Samples with presence of fibrin should be centrifuged. Do not use hemolyzed, lipemic or contaminated samples. In bibliography⁽¹⁾ it is reported the following stability in serum:
- Refrigerated: 7 days
 - Frozen: 90 days

PROCEDURE

To program and calibrate assays, follow the instructions for use of the analyzer used, with the recommended general parameters that are detailed below.

Detailed information on the parameters to program the assay on ADVIA® Chemistry XPT System and on ADVIA® 1800 Chemistry System is available in the documentation section of the website (www.3diag.com) or on request to the Customer Support Service (✉ support@3diag.com - ☎ +34 93 244 86 79).

Contact the Customer Support Service or the representative of Siemens Healthcare to program the assay on other analyzers of the ADVIA® Chemistry Systems series.

General Assay Parameters

- ①Dispense and mix:
 - * **BUF sTfR** 80 µl
 - * Sample/Control: 25 µl (diluted 1:20)
 - * Calibrator: 25 µl (diluted 1:5)
 - ②Incubate the mixture for the time set by the analyzer (approx. 5 minutes).
 - ③Dispense and mix:
 - * **REAG Ab sTfR** 20 µl
 - ④Incubate a fixed time between 30 and 60 seconds
 - ⑤Read absorbance A1 (Blank) at 596 nm
 - ⑥Incubate a fixed time of about 5 minutes
 - ⑦Read absorbance A2 (End Point) at 596 nm
 - ⑧Interpolate the absorbance increment (A2-A1) of the samples and controls in the curve obtained with the calibrators
- Samples with concentrations higher than the upper limit of the assay range should be analyzed again, by programming a larger sample dilution in the analyzer, or diluted manually with Physiological Solution, to recover a value close to the midpoint of the measurement range.

Calibration Parameters

- Use the **3diag - sTfR - CAL SET**.
- Discard the lower value calibrator (level 1) of the Set and replace it with Physiological Solution to perform the analyzer's Reagent Blank (zero concentration point).
- The values of the calibrators to be introduced in the analyzer must be the concentrations detailed in the table of assigned values of the Set multiplied by 4 (factor between the dilutions of the calibrators and the samples).
- It is recommended to program two replicates of each calibration point.

- The calibration is non-linear. For the calculation it is recommended to use a *Logit Log 3* adjustment.

The assay must be recalibrated, at least when a new batch of reagents is used or when its parameterization is changed.

PERFORMANCES OF THE METHOD

Detailed information on the characteristics and performances of the assay is given in the Technical Report, available on the website (www.3diag.com) or upon request to the Customer Support Service (✉ support@3diag.com - ☎ +34 93 244 86 79).

QUALITY CONTROL

To monitor performances, it is recommended that internal controls be inserted into each analytical series. It is recommended to use the controls of the **3diag - sTfR - CONTROL**.

In order to process the controls, it is necessary to deactivate, for calibrators and controls, the clot detection system.

Each laboratory should establish its own quality scheme and corrective actions if controls do not meet the assigned tolerances.

The reagents have been subjected to quality control checks and should react as described in these instructions. Therefore, as a general recommendation, in case the controls do not give the expected reaction, as a precaution all reagents should be considered unreliable until their operation has been checked.

TRACEABILITY

Values are referred to the *Reference Reagent Recombinant Soluble Transferrin Receptor (rsTfR)* (code: 07/202) of the *WHO (World Health Organization)*, and apply to free rsTfR monomer⁽²⁾.

REFERENCE INTERVALS

It is always advisable for each laboratory to establish its own reference values.

The bibliography⁽³⁾ reports reference values between 10.6 and 24.6 nmol/l, equivalent to 0.76 and 1.76 mg/l.

CLINICAL SIGNIFICANCE

sTfR is a good marker in the diagnosis of iron deficiency. sTfR enables evaluation of erythropoiesis without the need for cytological study of the bone marrow, as long as iron deficiency is excluded. Thus, it is helpful in managing response to erythropoietin treatment (EPO), and it has also been proposed for use in anti-doping control.

The concentration of sTfR increases significantly and at an early stage in patients suffering from iron deficiency. Because it does not act as an acute phase reactant, it is especially helpful in the differential diagnosis of iron-deficiency anaemia and secondary anaemia from chronic disorders (not caused by cancer) in the presence of acute or inflammatory conditions (commonly seen in elderly patients for example) which affect Ferritin measurements. Reports suggest that the sTfR/Log(Ft) ratio results in higher sensitivity and specificity of differentiation with respect to the individual magnitudes.

In its compilation of recommendations⁽⁴⁾, the *World Health Organization (WHO)* concludes that sTfR levels, combined with measurements of serum ferritin concentrations, provide the most effective assessment of a population's iron status.

SYMBOLS

In addition to the harmonized symbols provided on the European Standard EN 980:2008, in the labels and instructions of use has been used the complementary symbology proposed⁽⁵⁾ by the *EDMA (European Diagnostic Manufacturers Association)*, whose meaning is detailed below.

REAG	Reagent
Ab	Antibody / Antiserum
BUF	Buffer
sTfR	Soluble Transferrin Receptor

BIBLIOGRAPHY

- (1) *Mayo Medical Laboratories* (www.mayomedicallaboratories.com) website, date of consultation: 30th May 2017.
- (2) For further information visit *NIBSC (National Institute for Biological Standards and Control)* website: www.nibsc.org.
- (3) *Quest Diagnostics™* website (www.questdiagnostics.com), date of consultation: 17th June 2017.
- (4) "Serum transferrin receptor levels for the assessment of iron status and iron deficiency in populations. Vitamin and Mineral Nutrition Information System." Geneva: *World Health Organization*; 2014 (WHO/NMH/NHD/MNM/14.6; http://apps.who.int/iris/bitstream/10665/133707/1/WHO_NMH_NHD_EPG_14.6_eng.pdf?ua=1), date of consultation: 30th September 2019.
- (5) EDMA Labelling Task Force: "EDMA Symbols for IVD Reagents and Components - Revision, October 2009".

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
NOTE #1: Applications available for **Optilite®**, **Alinity c**, **Architect c**, **ADVIA®** series, **AU®** series, **cobas®** series, and other analyzers by request

IMAGE® and **AU®** are registered trademarks of Beckman Coulter, Inc, Fullerton, CA.
BN™ and **ADVIA®** are registered trademarks of Siemens Healthcare Diagnostics Products GmbH, Marburg, Deutschland.
Optilite® is a registered trademark of The Binding Site Group Ltd, Birmingham, U.K.
Alinity and **Architect** related brand marks are registered trademarks of Abbott Laboratories, Abbott Park, Illinois, USA
cobas® and related brand marks are trademarks of Roche Diagnostics Ltd, Rottkreuz, Switzerland

ALL PROTEINS INCLUDED IN  **Accuracy 365** Online Quality Control Further information and registration at www.accuracy365.com

Accuracy 365 aims to monitor and compare, between laboratories, the results obtained in the Internal Quality Control processes of end-users of specific protein control materials manufactured by **TRIMERO Diagnostics**.

Main features:

- 100 % free and lifetime.
- Cloud application, Available 24 hours, 365 days per year.
- Intuitive interface and very simple to use.
- Calculation and representation of statistical data in real time.
- Anonymous participation, data are only accessible by its owner.