

3diag - sTfR - 800

Soluble Transferrin Receptor - for *IMMAGE*[®] 800

REF TD-42681

IVD CE

ANNEX to Instruction for Use:

Application (NIPIA) for *IMMAGE*[®] and *IMMAGE*[®] 800 Immunochemical Systems

<u>PROTOCOL INFORMATION</u>			
Reagent Lot Number:	(^{*1})	Units:	mg/L or nmol/L
Cartridge Lot Number:	(^{*1})	Protocol:	Non-competitive NIPIA
Cartridge Serial Number:	(^{*1})	AGXS Limit:	not defined
Reagent Expiry date:	(^{*1})		
Sample or Dilution Volume:	10 µL	Diluent:	DIL1 (^{*2})
Reaction Buffer Volume:	195 µL	Buffer:	BUF sTfR (^{*3})
Calibration Dilution:	1 : 1 (neat)	Gain:	1
Compartment A Volume:	50 µL	Sample Dilution:	1 : 5 (^{*4})
Compartment B Volume:	0 µL (not used)	Reaction Time:	10.0 minutes
<u>CALIBRATION INFORMATION</u>			
Levels:	6	Replicates:	2
Update Level:	4 (^{*5})	Replicates:	2 (^{*5})
Model:	Logistic 4 parameters		
<u>NOTES</u>			
<p>(^{*1}): See the labels of the Reagent Cartridge. (^{*2}): Beckman Coulter P/N no. 447640. For the <i>IMMAGE</i>[®] it has to be defined as DIL10. (^{*3}): It has to be defined as BUF10 for the <i>IMMAGE</i>[®], and as BUF10 to BUF15 for the <i>IMMAGE</i>[®] 800. (^{*4}): It has to be programmed after performing the calibration process for the first time. (^{*5}): The use of the calibration update, with only one calibrator level, is discouraged.</p>			